

## Recommended protocol for OCT3/4

Obtained in run 53

20 Feb 2018

**Immunostainer**

Type: Dako Autostainer

**Primary antibody**

Clone: C10  
Producer: Santa Cruz  
Product no. / lot no.: SC-5279 / A2315  
Diluent: Da Vinci green  
Diluent producer / no.: Biocare / PD900M  
Dilution factor: 1:75  
Incubation time / temperature: 30 min. / 22°C

**Epitope retrieval, HIER**

Device: Microwave oven  
Buffer: Dako TRS High pH  
Heating time at max. temp.: 10 min.  
Maximum heating temp.: 100°C

**Visualization system**

Producer: Biocare  
Product / no: MACH4 / M4U534L  
Linker: MACH4 mouse probe  
Incubation time linker: 10 min.  
Incubation time polymer: 15 min.  
Incubation temperature: 22°C

**Chromogen**

Producer: Biocare  
Product / no: Betazoid DAB / BDB2004  
Incubation time / temperature: 5 min. / 22°C  
Enhancement: None

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.