

## Recommended protocol for CK-LMW

Obtained in run

20 Sep 2015

### Immunostainer

Type: Ventana Benchmark Ultra

### Primary antibody

Clone: CAM 5.2  
Producer: Becton Dickinson  
Product no. / lot no.: 345779 / 29678  
Format: Ready-To-Use (prediluted)  
Incubation time / temperature: 32 min. / 37°C

### Epitope retrieval, proteolysis

Enzyme: Protease 1  
Enzyme producer / no: Ventana / 760-2018  
Incubation time / temp: 8 min. / 37°C

### Visualization system

Producer: Ventana  
Product / no: UltraView Universal DAB Detection Kit / 760-500  
Incubation time polymer: 8 min.  
Incubation temperature: 36°C

### Chromogen

Producer: Ventana  
Product / no: ultraView Universa DAB Detection Kit / 760-500  
Incubation time / temperature: 8 min. / 36°C  
Enhancement: CuSO<sub>4</sub>

### Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.