

## Recommended protocol for CD68

Obtained in run 60

05 Aug 2020

**Immunostainer**

Type: Tissue Tek Genie

**Primary antibody**

Clone: GR021  
Producer: Genmed/Sakura Finetek  
Product no. / lot no.: 8275 / na  
Diluent: None  
Diluent producer / no.: None / None  
Dilution factor: 1:0  
Incubation time / temperature: 30 min. / 25°C

**Epitope retrieval, HIER**

Device: On Board / On Machine  
Buffer: Sakura Finetek Sakura-High pH Antigen retrieval buffer  
Heating time at max. temp.: 45 min.  
Maximum heating temp.: 98°C

**Visualization system**

Producer: Sakura Finetek  
Product / no: Tissue-Tek Genie Pro / 8826-K250  
Incubation time polymer: 20 min.  
Incubation temperature: 25°C

**Chromogen**

Producer: Sakura-Finetek  
Product / no: DAB / in detection kit 882  
Incubation time / temperature: 6 min. / 25°C  
Enhancement: CuSO4

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.