

## Recommended protocol for CD34

Obtained in run

26 Sep 2015

**Immunostainer**

Type: Dako Autostainer

**Primary antibody**

Clone: QBend10  
Producer: Dako  
Product no. / lot no.: M7165 / 00032597  
Diluent: Antibody Diluent  
Diluent producer / no.: Dako / S0809  
Dilution factor: 1:100  
Incubation time / temperature: 30 min. / 20°C

**Epitope retrieval, HIER**

Device: Microwave oven  
Buffer: Std Tris-EDTA / EGTA pH 9  
Heating time at max. temp.: 15 min.  
Maximum heating temp.: 100°C

**Visualization system**

Producer: Dako  
Product / no: EnVision FLEX / K8010/DM822  
Linker: None  
Incubation time polymer: 30 min.  
Incubation temperature: 20°C

**Chromogen**

Producer: Dako  
Product / no: DAB+ / DM827 / DM823 / K800  
Incubation time / temperature: 10 min. / 20°C  
Enhancement: None

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.