

## Recommended protocol for PAX8

Obtained in run 64

17 Feb 2022

**Immunostainer**

Type: Leica BOND Max

**Primary antibody**

Clone: BC12  
Producer: BioCare  
Product no. / lot no.: ACI438C / 050420A  
Diluent: Renoir Red Diluent  
Diluent producer / no.: Biocare Medical / PD904JJ  
Dilution factor: 1:25  
Incubation time / temperature: 25 min. / 24°C

**Epitope retrieval, HIER**

Device: On Board / On Machine  
Buffer: Leica Bond Epitope Retrieval Solution 2  
Heating time at max. temp.: 20 min.  
Maximum heating temp.: 99°C

**Visualization system**

Producer: Leica  
Product / no: Bond Refine / DS9800  
Incubation time linker: 8 min.  
Incubation time polymer: 8 min.  
Incubation temperature: 24°C

**Chromogen**

Producer: Leica  
Product / no: Bond Refine / DS9800  
Incubation time / temperature: 10 min. / 24°C  
Enhancement: None

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.