

## Recommended protocol for VIM

Obtained in run

27 Sep 2015

**Immunostainer**

Type: Autostainer, LabVision

**Primary antibody**

Clone: SP20  
Producer: NeoMarkers  
Product no. / lot no.: RM-9120-S / 9120S1001E  
Diluent: Antibody Diluent  
Diluent producer / no.: Dako / S3022  
Dilution factor: 1:400  
Incubation time / temperature: 30 min. / 20°C

**Epitope retrieval, HIER**

Device: Water bath  
Buffer: Dako TRS High pH  
Heating time at max. temp.: 20 min.  
Maximum heating temp.: 98°C

**Visualization system**

Producer: Dako  
Product / no: EnVision+ Dual Link System-HRP / K4061  
Incubation time polymer: 30 min.  
Incubation temperature: 20°C

**Chromogen**

Producer: Dako  
Product / no: DAB+ / K3467 / K3468  
Incubation time / temperature: 10 min. / 20°C  
Enhancement: None

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.