

## Recommended protocol for Napsin A

Obtained in run 66

22 Aug 2022

### Immunostainer

Type: Intellipath Biocare MEDical

### Primary antibody

Clone: BC15  
Producer: Biocare Medical  
Product no. / lot no.: API3043AA / 011222A  
Diluent: Da Vinci Green  
Diluent producer / no.: Biocare Medical / PD900L  
Dilution factor: 1:1  
Incubation time / temperature: 30 min. / 25°C

### Epitope retrieval, HIER

Device: PT-link / PT-module  
Buffer: Dako TRS High pH (3-1)  
Heating time at max. temp.: 20 min.  
Maximum heating temp.: 98°C

### Visualization system

Producer: Biocare MEDical  
Product / no: MACH-4UNIVERSAL HRP-POLYMER DETECTION SYSTEM  
Linker: MACH4 Mouse Probe  
Incubation time linker: 10 min.  
Incubation time polymer: 20 min.  
Amplifier: MACH4 HRP Polymer  
Incubation time amplifier: NA min.  
Incubation temperature: 25°C

### Chromogen

Producer: Biocare Medical  
Product / no: Betazoid DAB chromogen / BDB900F  
Incubation time / temperature: 5 min. / 25°C  
Enhancement: None

### Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided



## **Nordic Immunohistochemical Quality Control**

Department of Pathology, Aalborg University Hospital, Hospitalsbyen 5, DK-9260 Gistrup, Denmark

information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.