

## Recommended protocol for Napsin A

Obtained in run 66

29 Aug 2022

**Immunostainer**

Type: abcarta FAIP-48T

**Primary antibody**

Clone: 810B1C8  
Producer: abcarta  
Product no. / lot no.: PA102 / PA220520D07  
Diluent: Antibody Diluent  
Diluent producer / no.: abcarta / PB001  
Dilution factor: 1:0  
Incubation time / temperature: 15 min. / 30°C

**Epitope retrieval, HIER**

Device: On Board / On Machine  
Buffer: abcarta abcarta-EDTA/ER2 pH 9  
Heating time at max. temp.: 20 min.  
Maximum heating temp.: 100°C

**Visualization system**

Producer: abcarta  
Product / no: abcarta-HRP Polymer / PS300  
Linker: abcarta-Dual Linker  
Incubation time linker: 8 min.  
Incubation time polymer: 8 min.  
Incubation temperature: 30°C

**Chromogen**

Producer: abcarta  
Product / no: abcarta-DAB / PS300  
Incubation time / temperature: 5 min. / 30°C  
Enhancement: None

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.