

Recommended protocol for PAX8

Obtained in run 68

06 Feb 2023

Immunostainer

Type: Tissue-Tek Genie Advanced Stai

Primary antibody

Clone: RM436
Producer: Sakura Finetek
Product no. / lot no.: 8257-C010 / 8257A17922
Diluent: Sakura Finetek
Diluent producer / no.: Sakura Finetek / -
Dilution factor: 1:0
Incubation time / temperature: 30 min. / 20°C

Epitope retrieval, HIER

Device: Other
Buffer: Sakura Finetek Tissue-Tek Genie High pH Antigen Retrieval Buffer
Heating time at max. temp.: 45 min.
Maximum heating temp.: 98°C

Visualization system

Producer: Sakura Finetek
Product / no: Tissue-Tek Genie Pro Detection Kit DAB / 8826-K250
Linker: Tissue-Tek Genie Pro Detection Kit DAB
Incubation time linker: 15 min.
Incubation time polymer: 20 min.
Incubation temperature: 20°C

Chromogen

Producer: Sakura Finetek
Product / no: Tissue-Tek Genie Pro Detection Kit DAB / 8826-K250
Incubation time / temperature: 6 min. / 25°C
Enhancement: CuSO₄

Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.