

## Recommended protocol for URO II/III

Obtained in run 68

06 May 2023

**Immunostainer**

Type: Tissue-Tek Genie Advanced Stai

**Primary antibody**

Clone: BC21  
Producer: Sakura Finetek  
Product no. / lot no.: 8270-C010 / 8270A01234  
Diluent: Ab diluent  
Diluent producer / no.: Sakura Finetek / Sakura Finetek  
Dilution factor: 1:0  
Incubation time / temperature: 30 min. / 20°C

**Epitope retrieval, HIER**

Device: On Board / On Machine  
Buffer: Std Tris-EDTA / EGTA pH 9  
Heating time at max. temp.: 45 min.  
Maximum heating temp.: 98°C

**Visualization system**

Producer: Sakura Finetek  
Product / no: Tissue-Tek Genie Pro Detection Kit, DAB / 8826-K250  
Linker: Tissue-Tek Genie Pro Detection Kit, DAB / 8826-K250  
Incubation time linker: 15 min.  
Incubation time polymer: 20 min.  
Incubation temperature: 20°C

**Chromogen**

Producer: Sakura Finetek  
Product / no: Tissue-Tek Genie Pro Detection Kit, DAB / 8826-K250 / 8826-K250  
Incubation time / temperature: 6 min. / 25°C  
Enhancement: CuSO<sub>4</sub>

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.