

## Recommended protocol for Bcl-2

Obtained in run 73

05 Dec 2024

**Immunostainer**

Type: Celnovte CNT360

**Primary antibody**

Clone: C7B8  
Producer: Celnovte  
Product no. / lot no.: CBM-0041 / 202410001  
Diluent: None  
Diluent producer / no.: (None) / None  
Dilution factor: 1:0  
Incubation time / temperature: 20 min. / 25°C

**Epitope retrieval, HIER**

Device: On Board / On Machine  
Buffer: Std Tris-EDTA / EGTA pH 9  
Heating time at max. temp.: 20 min.  
Maximum heating temp.: 98°C

**Visualization system**

Producer: Celnovte  
Product / no: MicroStacker TM PLUS Polymer Detection Kit / SD5600  
Linker: MicroStacker TM PLUS Polymer Detection Kit  
Incubation time linker: 8 min.  
Incubation time polymer: 8 min.  
Incubation temperature: 25°C

**Chromogen**

Producer: Celnovte  
Product / no: MicroStacker TM PLUS Polymer Detection Kit / SD5600  
Incubation time / temperature: 6 min. / 25°C  
Enhancement: None

**Disclaimer:**

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.