**Nordic Immunohistochemical Quality Control**
Department of Pathology, Aalborg University Hospital, Hospitalsbyen 5, DK-9260 Gistrup, Denmark

# Recommended protocol for ER

Obtained in run B38

04 Jul 2024

**Immunostainer**

| | |
|---|---|
| Type: | Ventana Benchmark Ultra |

**Primary antibody**

| | |
|---|---|
| Clone: | SP1 |
| Producer: | Ventana Roche |
| Product no. / lot no.: | 790-4324/790-4325 / tbd |
| Format: | Ready-To-Use (prediluted) |
| Incubation time / temperature: | 16 min. / 36°C |

**Epitope retrieval, HIER**

| | |
|---|---|
| Device: | On Board / On Machine |
| Buffer: | Ventana CC1 |
| Heating time at max. temp.: | 64 min. |
| Maximum heating temp.: | 95°C |

**Visualization system**

| | |
|---|---|
| Producer: | Ventana |
| Product / no: | UltraView Universal DAB Detection Kit / 760-500 |
| Incubation time polymer: | 8 min. |
| Incubation temperature: | 36°C |

**Chromogen**

| | |
|---|---|
| Producer: | Ventana |
| Product / no: | ultraView Universal DAB Detection Kit / 760-500 |
| Incubation time / temperature: | 8 min. / 36°C |
| Enhancement: | CuSO4 |

**Disclaimer:**
NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.